

## A medium for rapid cleaning of *Aspergillus* cleistothecia

Susan G.W. Kaminskyj  
*Purdue University*

John E. Hamer  
*Purdue University*

Follow this and additional works at: <http://newprairiepress.org/fgr>

---

### Recommended Citation

Kaminskyj, S. G., and J.E. Hamer (1996) "A medium for rapid cleaning of *Aspergillus* cleistothecia," *Fungal Genetics Reports*: Vol. 43, Article 28. <https://doi.org/10.4148/1941-4765.1324>

This Short Note is brought to you for free and open access by New Prairie Press. It has been accepted for inclusion in *Fungal Genetics Reports* by an authorized administrator of New Prairie Press. For more information, please contact [cads@k-state.edu](mailto:cads@k-state.edu).

---

## A medium for rapid cleaning of *Aspergillus cleistothecia*

### **Abstract**

*Aspergillus cleistothecia* are usually cleaned of adherent Hülle cells by rolling them on an agar surface, which is a tedious process. We have found that incorporating 0.5 - 1.0 % diatomaceous earth (also called Celite) into 4% agar, as a mild abrasive, helps to speed cleaning about 5-fold (average 10 min reduced to 2 min).

### **Creative Commons License**



This work is licensed under a [Creative Commons Attribution-Share Alike 4.0 License](https://creativecommons.org/licenses/by-sa/4.0/).

## **A medium for rapid cleaning of *Aspergillus cleistothecia***

Susan G. W. Kaminskyj and John E. Hamer- Department of Biological Sciences, Hansen Research Building, Purdue University, West Lafayette, IN 47907-1392 USA

*Aspergillus cleistothecia* are usually cleaned of adherent Hülle cells by rolling them on an agar surface, which is a tedious process. We have found that incorporating 0.5 - 1.0 % diatomaceous earth (also called Celite) into 4% agar, as a mild abrasive, helps to speed cleaning about 5-fold (average 10 min reduced to 2 min). The medium may need to be stirred after autoclaving (stir bar sterilized separately to prevent boiling over) to ensure good mixing. Cleaning action is improved by pushing cleistothecium through a groove cut in the agar.

---

Last modified 7/25/96 KMC