

Method for replica-plating Neurospora

Marta Goodrich-Tanrikulu
United States Department of Agriculture

Thomas A. McKeon

Follow this and additional works at: <http://newprairiepress.org/fgr>

Recommended Citation

Goodrich-Tanrikulu, M., and T.A. McKeon (1994) "Method for replica-plating Neurospora," *Fungal Genetics Reports*: Vol. 41, Article 27. <https://doi.org/10.4148/1941-4765.1389>

This Regular Paper is brought to you for free and open access by New Prairie Press. It has been accepted for inclusion in Fungal Genetics Reports by an authorized administrator of New Prairie Press. For more information, please contact cads@k-state.edu.

Method for replica-plating Neurospora

Abstract

Method for replica-plating Neurospora

Creative Commons License



This work is licensed under a [Creative Commons Attribution-Share Alike 4.0 License](https://creativecommons.org/licenses/by-sa/4.0/).

Method for replica-plating *Neurospora*

Marta Goodrich-Tanrikulu and Thomas A. McKeon - United States Department of Agriculture, Agricultural Research Service, Western Regional Research Center, 800 Buchanan Street, Albany, CA 94710, U.S.A.

We have found a method to replica-plate colonies of *Neurospora crassa* which gives more reproducible results than velvet. We obtained exact replicas by overlaying plated conidia (on medium surface or in top agar of sorbose-containing media, i.e., modified Vogel's or synthetic crossing medium) with dry sterilized 0.45 m filters (Millipore HATF, Triton-free mixed cellulose acetate/nitrate). The colonies grow onto or through the filters, depending on the duration of contact. The filters can be transferred to new medium to make replicas, or used for other applications.