

Reciprocal translocation AR30 has a breakpoint distal to all known IIL markers.

J. F. Leslie

Follow this and additional works at: <https://newprairiepress.org/fgr>



This work is licensed under a [Creative Commons Attribution-Share Alike 4.0 License](https://creativecommons.org/licenses/by-sa/4.0/).

Recommended Citation

Leslie, J. F. (1982) "Reciprocal translocation AR30 has a breakpoint distal to all known IIL markers.," *Fungal Genetics Reports*: Vol. 29, Article 4. <https://doi.org/10.4148/1941-4765.1633>

This Research Note is brought to you for free and open access by New Prairie Press. It has been accepted for inclusion in Fungal Genetics Reports by an authorized administrator of New Prairie Press. For more information, please contact cads@k-state.edu.

Reciprocal translocation AR30 has a breakpoint distal to all known IIL markers.

Abstract

Reciprocal translocation AR30 has a breakpoint distal to all known IIL markers.

Reciprocal translocation AR30 has a breakpoint

distal to all known IIL markers.

T(IIL;VL)AR30 is a reciprocal translocation first described by Perkins and Barry (1977 Adv. Genet. 19: 113-285). Recently this rearrangement has been used in the synthesis of a complex chromosome rearrangement (Leslie, manuscript in preparation), and additional mapping data for both breakpoints were obtained. The II breakpoint is distal to all known markers on IIL, while the V breakpoint is distal to caf-1 on VL.

Three crosses were used to map the IIL breakpoint. Data from the cross T AR30 x Normal pi pyr-4 (Table I, cross A) place the AR30 breakpoint outside and a considerable distance from the pi pyr-4 region although the direction is not certain. The off-ratios in the first two classes can be attributed to decreased germination and/or growth of pi progeny.

TABLE 1
Crosses mapping AR30 breakpoint

				crossovers							
				Parentals	I	II	III	I&II	I&III	II&III	I, II & III
A.	I	II									
	N	pi	pyr-4								
	AR30	+	+	23	29	3		2			
	AR30	+	+	84	7	2		5			
B.	I										
	AR30	pyr-4	+	28	6						
	AR30	+	arg-5	30	12						
C.	I II III										
	AR30	+	ure-3 fl +	29	26	16	6	4	2	1	0
	AR30	arg-5	+	36	12	8	8	2	1	5	0
D.	I	II									
	N	caf-1	at	17	5	4		0			
	AR30	+	+	19	5	2		0			

- A. The number of normal (N) progeny for each class is given in the upper line.
 B. The number of AR30 progeny in the reciprocal classes is given in the lower line.
 C. The number of pyr-4 progeny for each class is given in the upper line.
 D. The number of arg-5 progeny in the reciprocal class is given in the lower line.
 E. The number of normal (N) progeny for each class is given in the upper line.
 F. The number of AR30 progeny in the reciprocal classes is given in the lower line.

Data from crosses homozygous for AR30 are also shown in Table I. In homozygous translocation crosses, markers distal to the IIL breakpoint will segregate independently of markers proximal to the IIL breakpoint. Markers on VL behave similarly. Cross B in Table I (AR30 pyr-4 x AR30 arg-5); indicates that pyr-4 and arg-5 remain linked so the AR30 breakpoint is not between pyr-4 and arg-5. Similarly, cross C shows that arg-5 is linked to ure-3, ure-3 is linked to fl, and fl is linked to rip: $P < 0.001$ from Chi-square tests for independence, for all four intervals. Thus, the II breakpoint of AR30 is on IIL distal to the known genetic markers but not terminal.

Data from the cross T AR30 x Normal at caf-1 Table I, cross D) show that the VL breakpoint is distal to caf-1. Perkins, Raju and Barry (1980 Chromosoma 76: 255-275) have shown cytologically that the AR30 breakpoint is proximal to the nucleolus organizer. Thus, the VL breakpoint of AR30 is between caf-1 and the nucleolus organizer.

A word of caution regarding this rearrangement is necessary. In heterozygous crosses (AR30 x normal sequence) many defective ascospores turn black. Instead of the 50% black spores expected from a reciprocal rearrangement crossed with normal sequence, about 85% of the shot spores from AR30 x Normal are black (the percentage of black spores increases with age). Consequently, the viability of the black spores from such a cross can be as low as 25 to 30%. In crosses homozygous for AR30 a normal situation prevails - 90 to 95% black spores are shot and 90 to 95% of the spores are viable. Thus, scoring AR30 visually by percent of shot black spores can be done only with practice, preferably by making comparisons of the tested stock when crossed with a known rearrangement tester and when crossed with a known normal-sequence tester. AR30 stocks containing fluffy are useful for such tests, and both T(IIL;VL)AR30 f1A (FGSC #3948 and T IIL:VL AR30 f1 a (FGSC #3949) have been deposited with FGSC. I have also deposited T(IIL;VL)AR30 caf-1 At a (FGSC #3950) and T(IIL;VL)AR30 caf-1 at a (FGSC #3951) which should be useful in detecting extremely distal IIL markers by linkage to caf-1 and at. (Supported by Public Health Service Grant AI-01462.) - - Department of Biological Sciences, Stanford University, Stanford, California 94305.

*Present address; Research and Development Division; International Minerals & Chemical Corporation, Terre Haute, Indiana 47808.

TABLE 1

Rates of growth and fungicide resistance of osmotic mutants of *Neurospora crassa*

Locus	Allele or isolation number	FGSC stock number	Rate of growth ^a				Resistance to fungicides ^b					
			MM	CM	MM + 2% NaCl	MM + 4% NaCl	procymidone	iprodione	vinclozolin	chloroneb	dicloran	quintozene
<u>os-1</u>	M16	812	50	19	0	5	+++	+	+++	+++	++	++
	M155-1	824	61	23	5	0	+++	+	+++	++	+	+
	B135	951	49	7	0	0	+++	++	+++	+++	++	+++
	P668	973	50	44	0	0	+++	+++	+++	+++	++	+++
	NM233(t)	1287	92	80	40	0	+++	+++	+++	++	++	+++
	P3282	1508	46	23	0	0	+++	+	+++	+++	++	+++
	NM204(t)	1637	80	60	9	0	+++	++	+++	++	++	++
	P3282	1644	49	11	0	0	+++	+	+++	+++	++	+++
	AR2	1675	63	21	0	0	+++	++	+++	++	++	+++
	P5990	2432	51	14	0	0	+++	++	+++	+++	++	+++
	P6549	2584	64	23	0	0	+++	++	+++	+++	++	+++
<u>os-2</u>	UCLA80	2238	106	109	62	8	++	++	++	+	++	+
	UCLA93	2240	98	95	45	0	++	++	++	+	++	+
<u>os-4</u>	NM2010	2429	106	112	58	32	++	++	++	+	+	+
<u>os-5</u>	NM2160	1638	98	98	47	0	++	++	++	+	+	+
<u>flm-1</u>	Y256M209	3624	94	76	48	0	+++	++	+++	++	++	+++
<u>flm-2</u>	Y256M223	2668	98	93	90	60	++	++	++	+	++	+
<u>cut</u>	LLM1	2385	96	102	43	12	0	0	0	0	0	0
<u>gla</u>	T9M150	3428	58	58	28	17	0	0	0	0	0	0
	T9M150	3429	108	116	92		0	0	0	0	0	0

^aIncrease in colony diameter, mm/24 hr, at 26°C; mean of 3 replicates.

CM = MM supplemented with 0.5% casamino acids and 0.5% yeast extract.

^b0 = sensitive (ED₉₅ < 2µg fungicide/ml); + = low resistance (ED₉₅ >2<10µg/ml); +++ = high resistance (ED₅₀ >50µg/ml; ED₉₅ >100µg/ml); ++ = intermediate levels of resistance. ED values are the concentrations of fungicide in agar media that reduce radial growth by 50% or 95%.